

BI+ BM300

PORTABLE MINI VEIN FINDER



www.portable-ultrasound-scanner.com

Quality, Professional, Innovative

BM300

PORTABLE MINI VEIN FINDER

Product Appearance



Product Features

- Six red light source, more powerful, no radiation, no laser
- Golden aluminum material shows high-end
- To show vascular in three-dimensional
- Brightness adjustment
- Easy to recharge with USB connection
- ergonomic design, lightweight, easy to carry, easy to operate

Product Applications

- Hospital, pediatric clinic, injection room
- Outpatient clinic, medical laboratory
- Blood bank, military medics unit, the CDC
- Rehabilitation center, nursing homes and others

Product Specifications

-Wavelength:	625nm
-Power in:	5V
-Brightness:	300lm
-Material:	golden aluminum

-Net weight:	0.15kg
-Length:	10cm
-Bottom diameter:	2.7cm

Wuxi Biomedical Technology Co., Ltd.

www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobilephone: 86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.